




Fitting Guide – Lower face cuff with link – fitting for a smaller patient

[YouTube Video – search Hilotherm Lower Face Cuff or scan QR code](#)

- Set up the device following the separate instructions
- Attach the cuff, turn the device on and wait for the cuff to inflate
- Select desired temperature (typically 16-19°C for facial treatment) and wait for a few minutes until the temperature has been achieved
- Lay out as shown below – Smooth side towards the patient's face



- Feed the longer hilofix strap through one of the slots at the top of the cuff, then down through the two slots at the very bottom of the cuff and up to the other slot at the top of the cuff to create a V shape.
- Attach the shorter strap to one of the slots in the middle of the cuff.
- The central cut out section should be in line with the mouth. Take the top ends of the longer hilofix strap and fix these on to each other on top of the head. The cut out section should fit around the mouth and the strap should sit under the chin as shown in the pictures.
- The smaller strap should run along the neck linking the 2 slots in the middle of the cuff
- Check the water is circulating all the way around the cuff by feeling the cuff by the input and output blue tubes



Hilotherm Ltd

Registered Office: Unit 3, The Venture Centre, University of Warwick Science Park,
Sir William Lyons Rd, Coventry CV4 7EZ

Office +44 (0) 2476 323 372 Fax +44 (0) 843 304 5071


Email: info@hilotherm.co.uk Registered in England: 11409331



Fitting Guide – Lower face cuff with link

YouTube Video – Search Hilotherm Lower Face Cuff or scan QR code

- Set up the device following the separate instructions
- Attach the cuff, turn the device on and wait for the cuff to inflate
- Select desired temperature (typically 16-19°C for facial treatment) and wait for a few minutes until the temperature has been achieved
- Lay out as shown below with the smooth side towards the patient's face
- The shorter hilofix bandage attached to the top right hand slot A and the other end through slot E to create a loop.



- Feed the second longer hilofix strap through slot B then C then D



- Place the loop strap (between slot A and E) on top of the head and line up the central cut out section with the mouth. Tighten the strap as needed to hold securely in place whilst keeping the open area around the mouth.
- Take either end of the second hilofix strap (by point B and D) and bring these to meet around the back of the neck. Cross the straps over and Velcro in place.
- Check the water is circulating all the way around the cuff by feeling the cuff by the input and output blue tubes

Hilotherm Ltd

Registered Office: Unit 3, The Venture Centre, University of Warwick Science Park,
Sir William Lyons Rd, Coventry CV4 7EZ

Office +44 (0) 2476 323 372 Fax +44 (0) 843 304 5071

Email: info@hilotherm.co.uk Registered in England: 11409331